



ABILITYER III

Microcatheter


***1.6 Fr Distal O.D. &
0.016" GW Compatibility***

HI-LEX CORPORATION

ABILITYER III


Super-selective approach enabled by 0.016" GW compatibility and 1.6 Fr distal tip!!

Peripheral deliverability





* The arrow indicates the microcatheter's distal tip.

Fluoroscopic visibility



* Imaged with an in-house X-ray CT system.

Flow rate



Maximum pressure resistance (distal end open): 1,000 psi (6,895 kPa)

Product No.	Effective length	Distal tip shape	O.D.			I.D.		Max. GW	DSV
			Distal tip	Distal	Proximal	Distal	Proximal		
HLAB3-16-105ST	105 cm	Straight	1.6 Fr	1.8 Fr (0.61 mm)	2.8 Fr (0.93 mm)	0.017" (0.43 mm)	0.023" (0.58 mm)	0.016"	0.40 mL
HLAB3-16-130ST	130 cm								0.46 mL
HLAB3-16-150ST	150 cm								0.51 mL
HLAB3-16-105AG	105 cm	45° Angled	1.6 Fr	1.8 Fr (0.61 mm)	2.8 Fr (0.93 mm)	0.017" (0.43 mm)	0.023" (0.58 mm)	0.016"	0.40 mL
HLAB3-16-130AG	130 cm								0.46 mL
HLAB3-16-150AG	150 cm								0.51 mL

Included Accessories : Shaping mandrel, Inserter

■ Medical Device: HL Microcatheter 3

■ EOG Sterilized, Disposable: Single use only

■ Japanese Marketing Approval Number: 30400BZX00041000

■ Please read the instructions for use prior to device application.

(All data presented in this catalog are derived from measurements conducted by our company.)

Manufacturer :

HI-LEX CORPORATION

Hi-Lex Corporation

1-12-28 Sakaemachi, Takarazuka, Hyogo 665-0845, Japan
 Tel: +81-797-85-2571 Fax: +81-797-85-2990